


## Mouse Anti-Rat IgG<sub>2c</sub> Clone 2C8F1 Biotin mAb

<b>Catalog No.</b>	CMR024	<b>Quantity:</b>	500 µg
<b>Description:</b>	Mouse Anti-Rat IgG <sub>2c</sub> Clone 2C8F1 Biotin mAb		
<b>Concentration:</b>	500 µg/mL		
<b>Specificity:</b>	Rat IgG <sub>2c</sub>		
<b>Immunogen:</b>	Rat IgG <sub>2c</sub> hybridoma		
<b>Isotype:</b>	Mouse (BALB/c) IgG <sub>2a</sub> K		
<b>Clone:</b>	2C8F1		
<b>Conjugate:</b>	Biotin		
<b>Formulation:</b>	Liquid volume of 1.0 mL in a buffer of PBS/NaN <sub>3</sub> .		
<b>Isotype Control</b>	Mouse IgG <sub>2a</sub> -BIOT		
<b>Applications:</b>	Applications for relevant formats of this clone include ELISA, FLISA, Flow Cytometry, and Immunohistochemistry.		
<b>Application Notes:</b>	Working dilution for ELISA 1:5,000 - 1:10,000, for Flow Cytometry ≤ 1 µg/10 <sup>6</sup> cells. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
<b>Storage &amp; Stability:</b>	Store at 2-8°C. For long term storage at -20°C, aliquot into smaller volumes to avoid multiple freeze/thaw cycles.		

ELISA plate was coated with purified rat IgG1, IgG2a, IgG2b, IgG2c, and IgM. Immunoglobulins were detected with serially diluted Mouse Anti-Rat IgG2c-BIOT followed by Streptavidin-HRP (Cat. No. CPO016).



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.