

Rat Anti-Mouse VEGFD mAb

Catalog No.	CMV107	Quantity:	100 µg
Description:	Rat Anti-Mouse VEGFD monoclonal antibody		
Alternate Names:	Figf, c-fos induced growth factor, VEGF-D		
Specificity:	Recognizes mouse VEGFD		
Host:	Rat		
Immunogen:	Mouse VEGFD recombinant protein		
Isotype:	IgG2		
Formulation:	Lyophilized from a 0.2 µm sterile filtered solution in PBS		
Purification:	Protein G affinity chromatography		
Reconstitution:	Centrifuge vial prior to opening. Reconstitute with sterile distilled water.		
Cross-Reactivity:	No cross-reactivity with mouse VEGFA and VEGFB		
Applications	Western Blot Immunohistochemistry		
	The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Store lyophilized antibody at -80°C. Reconstituted antibody is stable for at least six months in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

