

## Normal Cynomolgus Monkey Serum, Pooled Donors, Off the Clot

<b>Catalog No.</b>	CNK001A	<b>Quantity:</b>	50 mL
<b>Description:</b>	Normal Cynomolgus monkey serum is recovered from whole blood donations from normal, healthy cyno monkeys born in the USA. The serum is collected from spontaneously-clotted whole blood and centrifuged to separate the serum from the clot. Sera from $\geq 3$ donors are pooled for this product.		
<b>Volume:</b>	50 mL		
<b>Source:</b>	<i>Macaca fascicularis</i>		
<b>Storage &amp; Stability:</b>	Store in working aliquots at $-80^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.		
<b>Origin:</b>	USA		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

