

Goat Anti-Hamster IgG (H+L) Biotin pAb

Catalog No.	CPA004	Quantity:	1.0 mg
Description:	Goat Anti-Hamster IgG (H+L) Biotin pAb. Pooled antisera from goats hyperimmunized with hamster IgG. Purified by affinity chromatography on hamster IgG covalently linked to agarose.		
Concentration:	500 µg/mL		
Specificity:	Reacts with the heavy and light chains of hamster IgG; may only react with Syrian hamster IgG		
Immunogen:	Hamster IgG		
Isotype:	Goat IgG		
Conjugate:	HRP (Horseradish Peroxidase)		
Formulation:	Liquid volume of 2.0 mL supplied in PBS/NaN ₃ .		
Isotype Control	Goat IgG-Biotin		
Cross-Reactivity:	May react with immunoglobulins from other species and the light chains of other hamster immunoglobulins		
Applications:	Applications for this product include ELISA, FLISA, and Western Blot.		
Application Notes:	Working dilution for ELISA 1:5,000 - 1:20,000, for FLISA 1:200 - 1:400. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C. For long term storage at -20°C, aliquot into smaller volumes to avoid multiple freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

