

SERPINA1

Goat Anti-human Alpha-1 AntiTrypsin (A1AT) Affinity Purified pAb

| | | | |
|---------------------------------|---|------------------|----------------|
| Catalog No. | CPA104A CPA104B | Quantity: | 100 µg 1 mg |
| Alternate Names: | asPI; A1A; AAT; PI1; A1AT; PRO2275; alpha1AT | | |
| Description: | A1AT, Goat Anti-human A1AT Affinity Purified pAb A1AT is a secreted serine protease inhibitor whose targets include elastase, plasmin, thrombin, trypsin, chymotrypsin, and plasminogen activator. | | |
| Concentration: | Typically > 1.0 mg/ml | | |
| Specificity: | Human Alpha 1 Antitrypsin (AAT) | | |
| Host: | Goat | | |
| Immunogen: | Human Alpha 1 Antitrypsin (AAT) | | |
| Formulation: | Liquid | | |
| Storage & Stability: | Product is stable for 1 year when stored at 2-8°C. | | |

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

