

BCL3

Rabbit Anti-Human B-Cell CLL/Lymphoma 3 pAb

Catalog No.	CPB106	Quantity:	100 µg
Alternate Names:	B-cell leukemia/lymphoma 3, BCL-3		
Gene ID:	602		
Description:	Rabbit Anti-Human BCL3 polyclonal antibody. BCL3, a proto-oncogene, is believed to be involved in certain human B cell leukemias. It encodes a protein that functions as an IκB-like molecule for native NFκB, but is specific for the p50 subunit of NFκB.		
Concentration:	0.2 mg/ml		
Specificity:	Detects 55-60 kDa BCL3 of human, mouse, and rat origins		
Host:	Rabbit		
Immunogen:	Synthetic peptide surrounding amino acid 424 of human BCL3		
Formulation:	Liquid in PBS, pH 7.2 + 30% glycerol + 0.5% BSA + 0.01% thimerosal. Precaution: Thimerosal is a poisonous and hazardous substance which should be handled by trained staff only.		
Purification:	Protein G purified		
Cross-Reactivity:	Reacts with human, mouse, and rat.		
Applications:	Western Blot Immunoprecipitation Immunohistochemistry		
Application Notes:	For Western Blot, use a working dilution of 2-4 µg/ml. For Immunoprecipitation, use a working dilution of 10-20 µg/ml. For Immunohistochemistry, use a working dilution of 20-40 µg/ml. Blocking Peptide (Cat. No. CPB106BP) is available separately. The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Store at -20°C or in working aliquots at -80°C for long term storage. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
480 Neponset Street
Bldg 12A
Canton, MA 02021

Toll Free: 888-769-1246
Phone: 781-828-0610
Fax: 781-828-0542

E-mail: info@cellsciences.com
Website: www.cellsciences.com