

## S

### Human Anti-SARS-CoV-2 Spike-RBD (Clone REGN10933) Neutralizing mAb

<b>Catalog No.</b>	CPC511A CPC511B	<b>Quantity:</b>	50 µg 100 µg
<b>Alternate Names:</b>	Spike glycoprotein, S glycoprotein receptor binding domain, S-RBD		
<b>Description:</b>	<p>Recombinant Human anti-SARS-CoV-2 Spike Protein Receptor Binding Domain, Clone REGN10933, is expressed in XtenCHO.</p> <p>Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) is an enveloped, positive-sense, single-stranded RNA virus that causes coronavirus disease 2019 (COVID-19). The structural proteins of SARS-CoV-2 include the envelope protein (E), spike or surface glycoprotein (S), membrane protein (M) and the nucleocapsid protein (N). The spike glycoprotein is found on the outside of the virus particle and gives coronavirus viruses their crown-like appearance. Spike glycoprotein is cleaved into the following 3 chains, Spike protein S1, Spike protein S2, Spike protein S2'. Spike protein S1 attaches the virion to the cell membrane by interacting with host receptor, initiating the infection. Binding to human ACE2 receptor and internalization of the virus into the endosomes of the host cell induces conformational changes in the Spike glycoprotein. Surface glycoprotein is an important target for vaccine development, antibody therapies and diagnostic antigen-based tests.</p>		
<b>UniProt ID:</b>	P0DTC2		
<b>Immunogen:</b>	Recombinant SARS-CoV-2 Spike-RBD protein		
<b>Specificity:</b>	Recognizes SARS-CoV-2 Spike-RBD protein		
<b>Bioactivity:</b>	EC <sub>50</sub> = 66.32 ng/ml with SARS-CoV-2 Spike-RBD		
<b>Source:</b>	XtenCHO		
<b>Isotype:</b>	Human IgG		
<b>Clone:</b>	REGN10933		
<b>Concentration:</b>	1.0 mg/ml		
<b>Formulation:</b>	Sterile-filtered PBS, pH 7.5 preservative free.		
<b>Purification:</b>	Protein A affinity chromatography		
<b>Applications:</b>	Neutralizing ELISA: suggested dilution 1:5,000 - 1:10,000 Western blot: suggested dilution 1:1,000 - 1:2,000		
<b>Storage &amp; Stability:</b>	Stable at 2-8°C for 1 week or for up to 1 year at -20°C to -80°C. It is recommended to prepare working aliquots of undiluted product and store -20°C to -80°C.		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

