

Goat Anti-Chicken IgY(H+L)-HRP pAb

Catalog No.	CPF000	Quantity:	1.0 mL
Description:	Goat Anti-Chicken IgY(H+L)-HRP pAb. Pooled antisera from goats hyperimmunized with IgY. Purification by affinity chromatography on chicken IgY covalently linked to agarose.		
Specificity:	Reacts with the heavy and light chains of chicken IgY		
Immunogen:	Chicken IgY		
Isotype:	Goat IgG		
Conjugate:	HRP (Horseradish Peroxidase)		
Formulation:	Liquid volume of 1.0 mL in a stock solution in 50% glycerol/50% PBS, pH 7.4. No preservative added.		
Isotype Control	Goat IgG-HRP		
Cross-Reactivity:	No cross adsorption. May react with immunoglobulins from other species and the light chains of other chicken immunoglobulins		
Applications:	Applications for this product include ELISA, Immunocytochemistry, Western Blot and Surface Plasmon Resonance.		
Application Notes:	Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C or long-term at -20°C.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

