

Goat Anti-Chicken IgY(H+L)-Biotin pAb

Catalog No.	CPF002	Quantity:	1.0 mg
Description:	Goat Anti-Chicken IgY(H+L)-Biotin pAb. Pooled antisera from goats hyperimmunized with IgY. Purification by affinity chromatography on chicken IgY covalently linked to agarose.		
Concentration:	500 µg/mL		
Specificity:	Reacts with the heavy and light chains of chicken IgY		
Immunogen:	Chicken IgY		
Isotype:	Goat IgG		
Conjugate:	Biotin		
Formulation:	Liquid volume of 2.0 mL in a buffer of PBS/NaN ₃ .		
Isotype Control	Goat IgG-Biotin		
Cross-Reactivity:	No cross adsorption. May react with immunoglobulins from other species and the light chains of other chicken immunoglobulins		
Applications:	Applications for this product include ELISA, FLISA, Flow Cytometry, Immunocytochemistry, and Western Blot.		
Application Notes:	Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C. For long term storage at -20°C, aliquot into smaller volumes to avoid multiple freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

