

Sheep Anti-Human Factor VII pAb

Catalog No.	CPF120A CPF120B	Quantity:	1 mg 5 mg
Alternate Names:	Coagulation Factor VII, F7		
Description:	Sheep Anti-Human Factor VII polyclonal antibody		
Concentration:	5.0 mg/ml		
Specificity:	Human Factor VII		
Host:	Sheep		
Formulation:	Liquid in 50% glycerol		
Storage & Stability:	Store in working aliquots at -20°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

