

## Goat Anti-Human Ig Fc (Ads) HRP pAb

<b>Catalog No.</b>	CPH016	<b>Quantity:</b>	1.0 mL
<b>Description:</b>	Goat Anti-Human Ig Fc (Ads) HRP. Pooled antisera from goats hyperimmunized with human IgG, IgM and IgA. Cross absorbed on human IgG Fab. Purified by affinity chromatography.		
<b>Specificity:</b>	Reacts with the Fc region of human IgG		
<b>Immunogen:</b>	Human IgG, IgM, IgA		
<b>Isotype:</b>	Goat IgG		
<b>Conjugate:</b>	HRP (Horseradish Peroxidase)		
<b>Formulation:</b>	Liquid volume of 1.0 mL supplied in a stock solution in 50% glycerol/50% PBS, pH 7.4. No preservative added.		
<b>Isotype Control</b>	Goat IgG-HRP		
<b>Cross-Reactivity:</b>	May react with IgG from other species.		
<b>Applications:</b>	Applications for this product include ELISA.		
<b>Application Notes:</b>	Working dilution for ELISA 1:4,000 - 1:8,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
<b>Storage &amp; Stability:</b>	Store at 2-8°C or long-term at -20°C.		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**