


Goat anti-Human IgG (IgA, IgM ads.) Affinity Purified pAb

Catalog No.	CPH028	Quantity:	1.0 mg
Description:	Pooled antisera from goats hyperimmunized with human IgG. Cross absorbed on human IgA and IgM. Purified by affinity chromatography using human IgG covalently linked to agarose.		
Specificity:	Reacts with the heavy chain of human IgG		
Immunogen:	Human IgG		
Isotype:	Goat IgG		
Purification:	Cross absorbed on human IgA and IgM. Affinity chromatography on human IgG covalently linked to agarose		
Concentration:	1.0 mg/ml purified IgG		
Formulation:	Borate buffered saline, pH 8.2 without carrier or preservatives.		
Cross-Reactivity:	May react with IgG from other species.		
Applications:	ELISA, Western blot		
Storage & Stability:	Store unopened at 2-8°C for up to 1 year.		

ELISA plate was coated with purified human IgG, IgM, and IgA. Immunoglobulins were detected with serially diluted Goat Anti-Human IgG (CPH028) followed by Mouse Anti-Goat IgG Fc-HRP



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

