

## Goat Anti-Rhesus IgG (H+L)-HRP

<b>Catalog No.</b>	CPK002	<b>Quantity:</b>	1.0 mL
<b>Description:</b>	Goat Anti-Rhesus IgG (H+L)-HRP. Pooled antisera from goats hyperimmunized with rhesus IgG. Purified by affinity chromatography on rhesus IgG covalently linked to agarose.		
<b>Specificity:</b>	Reacts with the heavy and light chains of rhesus IgG		
<b>Isotype:</b>	Goat IgG		
<b>Conjugate:</b>	HRP (Horseradish Peroxidase)		
<b>Formulation:</b>	Liquid volume of 1.0 mL in a stock solution of 50% glycerol, 50% PBS, pH 7.4. No preservative added.		
<b>Cross-Reactivity:</b>	May react with immunoglobulins from other species and the light chains of other rhesus immunoglobulins		
<b>Applications:</b>	ELISA, Immunocytochemistry, Western Blot		
<b>Application Notes:</b>	Working dilution for ELISA 1:4,000 - 1:8,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
<b>Storage &amp; Stability:</b>	Store at 2-8 °C or long term at -20 °C		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

