

## F12

### Rabbit Anti-Mouse Factor XII Affinity Purified pAb

<b>Catalog No.</b>	CPM008A CPM008B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	HAF, Hageman factor, FXII		
<b>Description:</b>	Rabbit affinity purified polyclonal IgG fraction to mouse factor XII.		
<b>Gene ID:</b>	58992		
<b>Specificity:</b>	Recognizes mouse Factor XII		
<b>Host:</b>	Rabbit		
<b>Immunogen:</b>	Recombinant full length wild-type protein from insect cells (glycosylated)		
<b>Isotype:</b>	IgG		
<b>Formulation:</b>	Frozen liquid in 0.05M Sodium Phosphate + 0.1 M NaCl + 1 mM EDTA, pH 6.6.		
<b>Purification:</b>	Affinity purified by immobilized mouse factor XII		
<b>Applications:</b>	ELISA, WB		
<b>Application Notes:</b>	The optimal concentration should be determined by the user for each specific application.		
<b>Storage &amp; Stability:</b>	When stored at -80°C, product is stable for 5 years from date of delivery. <b>Avoid repeated freeze-thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**