


Goat Anti-Mouse Kappa (Ads) HRP pAb

Catalog No.	CPM023	Quantity:	1.0 mL
Description:	Goat Anti-Mouse Kappa (Ads) HRP pAb. Pooled antisera from goats hyperimmunized with mouse κ light chains. Cross adsorbed on mouse λ light chains. Purified by affinity chromatography on mouse κ light chains covalently linked to agarose.		
Specificity:	Reacts with mouse κ light chains		
Immunogen:	Mouse κ light chains		
Isotype:	Goat IgG		
Conjugate:	HRP (Horseradish Peroxidase)		
Formulation:	Liquid volume of 1.0 mL supplied in a stock solution in 50% glycerol/50% PBS, pH 7.4. No preservative added.		
Isotype Control	Goat IgG-HRP		
Cross-Reactivity:	May react with κ light chains from other species		
Applications:	Applications include ELISA, ELISPOT, Immunohistochemistry, Immunocytochemistry, and Western Blot. .		
Application Notes:	Working dilution for ELISA 1:2,000 - 1:8,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C or long-term at -20°C.		

Mouse immunoglobulins were resolved by electrophoresis under reducing conditions, transferred to PVDF membrane, and probed with Goat Anti-Mouse Kappa-HRP. Proteins were visualized using chemiluminescent detection.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

