


Goat Anti-Mouse IgG_{2b} (Ads) HRP pAb

Catalog No.	CPM033	Quantity:	1.0 mL
Description:	Goat Anti-Mouse IgG _{2b} (Ads) HRP pAb. Pooled antisera from goats hyperimmunized with mouse IgG _{2b} . Cross adsorbed on mouse IgG ₁ , IgG _{2a} , IgG ₃ , IgM, and IgA; human immunoglobulins and pooled human sera. Purified by affinity chromatography on mouse IgG _{2b} covalently linked to agarose.		
Specificity:	Reacts with the heavy chain of mouse IgG _{2b}		
Immunogen:	Mouse IgG _{2b}		
Isotype:	Goat IgG		
Conjugate:	HRP (Horseradish Peroxidase)		
Formulation:	Liquid volume of 1.0 mL supplied in a stock solution in 50% glycerol/50% PBS, pH 7.4. No preservative added.		
Isotype Control	Goat IgG-HRP		
Cross-Reactivity:	May react with immunoglobulins from other species		
Applications:	Applications include ELISA, ELISPOT, Immunohistochemistry, Immunocytochemistry, and Western Blot.		
Application Notes:	Working dilution for ELISA 1:4,000 - 1:8,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C or long-term at -20°C.		

ELISA plate was coated with purified mouse IgG1, IgG2a, IgG2b, IgG3, IgM, and IgA. Immunoglobulins were detected with serially diluted Goat Anti-Mouse IgG2b HRP.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

