

Sheep Anti-FITC-Biotin

Catalog No.	CPO011	Quantity:	1.0 mg
Description:	Sheep Anti-FITC-Biotin. Pooled antisera from sheep hyperimmunized with fluorescein isothiocyanate (Isomer I).		
Concentration:	500 µg/mL		
Specificity:	Reacts with Fluorescein Isothiocyanate (FITC)		
Isotype:	Sheep IgG		
Conjugate:	Biotin		
Formulation:	Liquid volume 2.0 mL in a buffer of PBS/NaN ₃ .		
Isotype Control:	Sheep IgG-Biotin		
Cross-Adsorption:	None; may react with other fluorescein derivatives		
Applications:	ELISA		
Application Notes:	Working dilution; ELISA 1:5,000 - 1:20,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C. For long term storage at -20°C, aliquot into smaller volumes to avoid multiple freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.