

Goat Anti-Porcine IgG (H+L) HRP pAb

Catalog No.	CPP001	Quantity:	1.0 mL
Description:	Goat Anti-Porcine IgG (H+L) HRP pAb. Pooled antisera from goats hyperimmunized with porcine IgG. Purified by affinity chromatography on porcine IgG covalently linked to agarose.		
Specificity:	Reacts with the heavy and light chains of porcine IgG		
Immunogen:	Porcine IgG		
Isotype:	Goat IgG		
Conjugate:	HRP (Horseradish Peroxidase)		
Formulation:	Liquid volume of 1.0 mL supplied in a stock solution in 50% glycerol/50% PBS, pH 7.4. No preservative added.		
Isotype Control	Goat IgG-HRP		
Cross-Reactivity:	May react with immunoglobulins from other species and the light chains of other porcine immunoglobulins.		
Applications:	Applications include ELISA, Immunocytochemistry, immunohistochemistry and Western Blot.		
Application Notes:	Working dilution for ELISA 1:4,000 - 1:8,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C or long-term at -20°C.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

