

Rabbit Anti-Human PAI-1 Affinity Purified pAb

Catalog No.	CPP107B CPP107C	Quantity:	100 µg 1 mg
Alternate Names:	Plasminogen Activator Inhibitor type 1, SERPINE1, PAI, PAI1, PLANH1		
Description:	Rabbit Anti-Human Plasminogen Activator Inhibitor type 1 Affinity Purified polyclonal antibody.		
Concentration:	10.0 mg/ml		
Specificity:	Human PAI-1		
Host:	Rabbit		
Isotype:	IgG		
Formulation:	Liquid in 0.05 M sodium phosphate + 0.1 M NaCl + 1 mM EDTA, pH 6.6		
Purification:	Biospecific affinity chromatography		
Applications:	Western Blot ELISA Immunohistochemistry The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Store in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

