

Goat Anti-Rabbit PAI-1 Affinity Purified pAb

| | | | |
|---------------------------------|--|------------------|----------------|
| Catalog No. | CPP125B CPP125C | Quantity: | 100 µg 1 mg |
| Alternate Names: | Plasminogen Activator Inhibitor type 1, SERPINE1, PAI, PAI1, PLANH1 | | |
| Description: | Goat Anti-Rabbit Plasminogen Activator Inhibitor type 1 Affinity Purified polyclonal antibody. | | |
| Concentration: | 10.0 mg/ml | | |
| Specificity: | Rabbit PAI-1 | | |
| Host: | Goat | | |
| Isotype: | IgG | | |
| Formulation: | Liquid in 0.05 M sodium phosphate + 0.1 M NaCl + 1 mM EDTA, pH 6.6 | | |
| Purification: | Biospecific affinity chromatography | | |
| Applications: | Western Blot ELISA Immunohistochemistry The optimal concentration should be determined by the user for each specific application. | | |
| Storage & Stability: | Store in working aliquots at -80°C. Avoid repeated freeze-thaw cycles. | | |

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

