

Donkey Anti-Rat IgG (H+L) (Ads) HRP pAb

Catalog No.	CPR012	Quantity:	1.0 mL
Description:	Donkey Anti-Rat IgG (H+L) (Ads) HRP pAb. Pooled antisera from donkeys hyperimmunized with rat IgG. Cross adsorbed on mouse serum proteins. Purified by affinity chromatography on rat IgG covalently linked to agarose.		
Specificity:	Reacts with the heavy and light chains of rat IgG and the light chains of rat IgM		
Immunogen:	Rat IgG		
Isotype:	Donkey IgG		
Conjugate:	HRP (Horseradish Peroxidase)		
Formulation:	Liquid volume of 1.0 mL in a stock solution of 50% glycerol, 50% PBS, pH 7.4. No preservative added.		
Isotype Control:	Donkey IgG-HRP		
Cross-Reactivity:	May react with immunoglobulins from other species and the light chains of other rat immunoglobulins		
Applications:	ELISA, Immunohistochemistry, Western Blot		
Application Notes:	Working dilution for ELISA 1:4,000 - 1:8,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8 °C or long term at -20 °C		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

