

## Sheep Anti-Human tPA Antiserum

<b>Catalog No.</b>	CPT127B CPT127C	<b>Quantity:</b>	1 ml 10 ml
<b>Alternate Names:</b>	Tissue Plasminogen Activator, T-Plasminogen Activator, PLAT, T-PA		
<b>Description:</b>	Sheep Anti-Human Tissue Plasminogen Activator Antiserum.		
<b>Concentration:</b>	Whole Antiserum		
<b>Specificity:</b>	Human tPA		
<b>Host:</b>	Sheep		
<b>Storage &amp; Stability:</b>	Store in working aliquots at -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

