

Rabbit Anti-Mouse tPA Biotin pAb

Catalog No.	CPT131B CPT131C	Quantity:	1 mg 10 mg
Alternate Names:	Tissue Plasminogen Activator, T-Plasminogen Activator, PLAT, T-PA		
Description:	Rabbit Anti-Mouse Tissue Plasminogen Activator Biotin labeled polyclonal antibody.		
Concentration:	10 mg/ml		
Specificity:	Mouse tPA		
Host:	Rabbit		
Isotype:	IgG		
Conjugate:	Biotin		
Formulation:	Liquid in 0.05 M sodium phosphate + 0.1 M NaCl + 1 mM EDTA, pH 6.6		
Purification:	Protein A purification		
Applications:	Western Blot ELISA The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Store in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

