

Sheep Anti-Multispecies tPA Antiserum

Catalog No.	CPT134B CPT134C	Quantity:	1 ml 10 ml
Alternate Names:	Tissue Plasminogen Activator, T-Plasminogen Activator, PLAT, T-PA		
Description:	Sheep Anti-Multispecies Tissue Plasminogen Activator Antiserum.		
Concentration:	Whole Antiserum		
Specificity:	Mouse and Rat tPA		
Host:	Sheep		
Storage & Stability:	Store in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

