

Rabbit Anti-Human Urokinase pAb

Catalog No.	CPU104B CPU104C	Quantity:	1 mg 10 mg
Alternate Names:	Urokinase Plasminogen Activator, PLAU, ATF, UPA, URK, u-PA		
Description:	Rabbit Anti-Human Urokinase polyclonal antibody		
Concentration:	10 mg/ml		
Specificity:	Human urokinase		
Host:	Rabbit		
Isotype:	IgG		
Formulation:	Liquid in 0.05 M sodium phosphate + 0.1 M NaCl + 1 mM EDTA, pH 6.6		
Purification:	Protein A purification		
Applications:	Western Blot ELISA Immunohistochemistry The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Store in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

