

Rabbit Anti-Mouse Urokinase Antiserum

| | | | |
|---------------------------------|--|------------------|---------------|
| Catalog No. | CPU111B CPU111C | Quantity: | 1 ml 10 ml |
| Alternate Names: | Urokinase Plasminogen Activator, PLAU, ATF, UPA, URK, u-PA | | |
| Description: | Rabbit Anti-Mouse Urokinase Antiserum | | |
| Concentration: | Whole Antiserum | | |
| Specificity: | Mouse Urokinase | | |
| Host: | Rabbit | | |
| Formulation: | Liquid | | |
| Applications: | Western Blot ELISA Immunohistochemistry The optimal concentration should be determined by the user for each specific application. | | |
| Storage & Stability: | Store in working aliquots at -80°C. Avoid repeated freeze-thaw cycles. | | |

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

