

## VTN

### Rabbit Anti-Human Vitronectin pAb

<b>Catalog No.</b>	CPV121B CPV121C	<b>Quantity:</b>	1 mg 10 mg
<b>Alternate Names:</b>	V75, VN, VNT		
<b>Gene ID:</b>	7448		
<b>Description:</b>	Rabbit Anti-Human Vitronectin polyclonal antibody		
<b>Concentration:</b>	10 mg/ml		
<b>Specificity:</b>	Human vitronectin		
<b>Host:</b>	Rabbit		
<b>Isotype:</b>	IgG		
<b>Formulation:</b>	Liquid in 0.05 M sodium phosphate + 0.1 M NaCl + 1 mM EDTA, pH 6.6		
<b>Purification:</b>	Protein A purification		
<b>Applications:</b>	Western Blot ELISA Immunohistochemistry The optimal concentration should be determined by the user for each specific application.		
<b>Storage &amp; Stability:</b>	Store in working aliquots at -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

