

Recombinant Human CK2 alpha 1 polypeptide GST Active

Catalog No.	CRC110A CRC110B	Quantity:	5 µg 10 µg
Alternate Names:	Casein Kinase 2 alpha 1 polypeptide, CSNK2A1, ACK2, CK II, CK2A1, CK2-alpha, CKII, CKII alpha		
Description:	Recombinant full length human CK2 alpha 1 polypeptide containing N-terminal GST tag.		
Concentration:	0.1 mg/ml		
GeneID:	1457		
Protein Accession No:	NM_001895		
Source:	Sf9 insect cells		
Molecular Weight:	~ 70 kDa		
Formulation:	Liquid in 50 mM Tris-HCl, pH 7.5 + 150 mM NaCl + 0.25 mM DTT + 0.1 mM EGTA + 0.1 mM EDTA + 0.1 mM PMSF + 25% glycerol		
Purity:	> 95% as determined by densitometry		
Storage & Stability:	Stable for 1 year in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.