

Native Human D-Dimer

Catalog No.	CRD111A CRD111B	Quantity:	100 µg 1.0 mg
Alternate Names:	Fragment D-dimer, Fibrin degradation product		
Description:	D-dimer is a fibrin degradation product (FDP), a small protein fragment present in the blood after a blood clot is degraded by fibrinolysis. It is so named because it is comprised of two D fragments of the fibrin protein joined by a cross-link. The concentration of D-dimer may be determined by a blood test in order to help diagnose thrombosis.		
Source:	Human plasma		
Appearance:	Clear colorless solution		
Concentration:	1-5 mg/ml, by $E^{0.1\%}_{280nm} = 1.78$		
Molecular Weight:	~200 kDa		
Formulation:	Sterile filtered solution of 0.05 M sodium phosphate, 0.15 M sodium chloride, 0.09% sodium azide, pH 7.5		
Purity:	≥95% by SDS-PAGE		
Identity:	Confirmed by ELISA		
Storage & Stability:	Centrifuge vial briefly to recover the entire contents. Store in working aliquots at -20°C to -80°C for up to 1 year. Avoid repeated freeze-thaw cycles.		
Statement:	Non-Infectious Disease Certification: Non-reactive for HIV-1/HCV/HBV by NAT; HBsAg, HCV Ab, HIV-1&2 Ab and RPR by currently approved FDA methods. However, because no test method can offer complete assurance that infectious agents are absent, this material should be handled at Bio-Safety Level 2 (BSL 2) as recommended for potentially infectious human serum or blood specimen in the CCD/NIH manual "Biosafety in Microbiological and Biomedical Laboratories", 2009.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com