

## Recombinant Human Eotaxin

<b>Catalog No:</b>	CRE000C	<b>Size:</b>	1.0 mg
<b>Lot Number:</b>	L122202		
<b>Molecular Weight:</b>	8.3 kDa		
<b>Purity:</b>	>98% as determined by SDS-PAGE analysis.		
<b>Biological Activity:</b>	ED <sub>50</sub> range = 5-20 ng/mL, determined by the dose dependent chemotaxis of human peripheral eosinophils. The optimal concentration should be determined for each specific application. Eotaxin is a $\beta$ -chemokine (C-C chemokine).		
<b>Formulation:</b>	Lyophilized, carrier free.		
<b>Sterility:</b>	Filtered through a 0.22 micron filter prior to lyophilization.		
<b>Endotoxin:</b>	<0.1 ng/ $\mu$ g		
<b>Source:</b>	Recombinant human eotaxin is produced in <i>E. coli</i> and purified via sequential chromatography.		
<b>Reconstitution:</b>	We recommend that the vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute lyophilized recombinant human eotaxin in sterile, distilled water to a concentration of 0.1-1.0 mg/mL. Further dilutions should be made in low endotoxin medium or a buffered solution containing a carrier protein such as heat inactivated FCS or tissue culture grade BSA.		
<b>Suggested Working Dilutions:</b>	The optimal concentration should be determined for each specific application.		
<b>Storage:</b>	Store lyophilized recombinant human eotaxin at 2-8°C, preferably desiccated. Upon reconstitution, apportion into working aliquots and store at $\leq$ -20°C. Avoid repeated freeze/thaw cycles.		
<b>References:</b>	<p>Ponath, P.D., et al. (1996) Cloning of the human eosinophil chemoattractant, eotaxin. Expression, receptor binding, and functional properties suggest a mechanism for the selective recruitment of eosinophils. <i>J. Clin. Invest.</i> 97:604-612.</p> <p>Garcia-Zepeda, E.A., et al. (1996) Human eotaxin is a specific chemoattractant for eosinophil cells and provides a new mechanism to explain eosinophilia. <i>Nat. Med.</i> 2:449-456.</p> <p>Kampen, G.T., et al. (2000) Eotaxin induces degranulation and chemotaxis of eosinophils through the activation of ERK2 and p38 mitogen-activated protein kinases. <i>Blood</i> 95:1911-1917.</p> <p>Van Drenth, C., et al. (2000) Desensitization of CXC chemokine receptor 4, mediated by IL-16/CD4, is independent of p56 (lck) enzymatic activity. <i>J. Immunol.</i> 165 (11):6356-6363.</p>		

**FOR RESEARCH USE ONLY. NOT FOR USE IN HUMAN DIAGNOSTIC OR THERAPEUTIC PROCEDURES.**



**Cell Sciences, Inc.**  
480 Neponset Street  
Building 12A  
Canton, MA 02021

Toll Free: 888 769-1246  
Phone: 781 828-0610  
Fax: 781 828-0542

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Web Site: [www.cellsciences.com](http://www.cellsciences.com)