

EFNA4

Recombinant Human Ephrin A4 GST

Catalog No.	CRE111A CRE111B CRE111C	Quantity:	2 µg 5 µg 10 µg
Alternate Names:	Ephrin-A4, EPH-related receptor tyrosine kinase ligand 4, LERK-4, EFNA4, EPLG4, LERK4, EFL4		
Description:	<p>Ephrin A4 is a member of the ephrin (EPH) family. Ephrins and EPH-related receptors comprise the largest subfamily of receptor protein-tyrosine kinases. Based upon their structures and sequence relationships, ephrins are divided into the ephrin-A (EFNA) class, which are attached to the membrane by a glycosylphosphatidylinositol linkage, and the ephrin-B (EFNB) class, which are transmembrane proteins. Class A ephrins are linked to the membrane by a GPI linkage and bind primarily to EphA receptors; Class B ephrins contain a membrane-spanning region and bind primarily to EphB receptors. Both ephrins and Eph receptors are largely expressed throughout the ectoderm, mesoderm, and endoderm of vertebrate embryos.</p> <p>EFNA binds to the EPHA2, EPHA4, EPHA5, EPHA6, and EPHA7 receptors. Two transcript variations which encode different isoforms were identified through sequence analysis. During development the Eph receptor and their membrane-attached ligands, ephrins, show diverse expression patterns.</p> <p>The ephrins and EPH-related receptors mediate developmental events, principally in the nervous system and in erythropoiesis. The critical roles of Eph and ephrins during development imply involvement of these genes in congenital disorders affecting the nervous system and other tissues.</p> <p>Ephrins and Eph receptors are possibly involved in carcinogenesis.</p> <p>Ephrin A4 Human Recombinant full length protein expressed in <i>E.coli</i>, shows a 46 kDa SDS-PAGE (Including GST tag). The Ephrin A4 is purified by proprietary chromatographic techniques.</p>		
Concentration:	0.1 mg/ml		
Gene ID:	1945		
Source:	<i>E. coli</i>		
Molecular Weight:	46 kDa		
Formulation:	Sterile filtered liquid in 50 mM Tris-Acetate, pH 7.5 + 1 mM EDTA + 20% glycerol		
Applications:	<p>Western Blot ELISA Inhibition Assays</p> <p>The optimal concentration should be determined by the user for each specific application.</p>		
Storage & Stability:	Stable for 1 year at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



® Cell Sciences®
480 Neponset Street
Bldg 12A
Canton, MA 02021

Toll Free: 888-769-1246
Phone: 781-828-0610
Fax: 781-828-0542

E-mail: info@cellsciences.com
Website: www.cellsciences.com