

## Native Bovine Fibrinogen

<b>Catalog No.</b>	CRF138B CRF138C	<b>Quantity:</b>	25 mg 100 mg
<b>Description:</b>	Native Bovine Fibrinogen		
<b>Concentration:</b>	40 mg/ml		
<b>Source:</b>	Bovine plasma		
<b>Molecular Weight:</b>	330 kDa		
<b>Formulation:</b>	Liquid in 0.02 M sodium citrate-HCl, pH 7.4		
<b>Storage &amp; Stability:</b>	Store in working aliquots at -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

