

Native Human Factor Xa beta

Catalog No.	CRF147B CRF147C	Quantity:	500 µg 1 mg
Description:	Native Human Factor Xa beta		
Concentration:	5.0 mg/ml		
Source:	Human plasma		
Molecular Weight:	46 kDa		
Formulation:	Liquid in 50% (v/v) glycerol/H ₂ O		
Specific Activity:	Measured in a Factor Xa clotting assay and/or chromogenic substrate assay.		
Storage & Stability:	Store in working aliquots at -20°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

