


TNFRSF11B

Recombinant Human Osteoprotegerin / TNFRSF11B (Fc Tag)

Catalog No.	CRH435A-Fc CRH435B-Fc	Quantity:	50 µg 100 µg
Alternate Names:	Tumor necrosis factor receptor superfamily member 11B, Osteoclastogenesis inhibitory factor, Osteoprotegerin		
Description:	<p>Osteoprotegerin or TNFRSF11B is a member of the TNF-receptor superfamily. This protein is an osteoblast-secreted decoy receptor that functions as a negative regulator of bone resorption. This protein specifically binds to its ligand, osteoprotegerin ligand, both of which are key extracellular regulators of osteoclast development. Studies of the mouse counterpart also suggest that this protein and its ligand play a role in lymph-node organogenesis and vascular calcification. Alternatively spliced transcript variants of this gene have been reported, but their full length nature has not been determined.</p> <p>Osteoprotegerin/TNFRSF11B acts as decoy receptor for RANKL and thereby neutralizes its function in osteoclastogenesis. This protein may inhibit the activation of osteoclasts and promotes osteoclast apoptosis in vitro. Bone homeostasis seems to depend on the local RANKL/OPG ratio. Osteoprotegerin/TNFRSF11B also play a role in preventing arterial calcification, act as decoy receptor for TRAIL and protect against apoptosis. TRAIL binding blocks the inhibition of osteoclastogenesis.</p>		
UniProt ID:	O00300		
Accession Number:	NP_002537.3		
Protein Construction:	The DNA sequence encoding the human TNFRSF11B (Met 1- Leu 401) was expressed with the Fc region of human IgG1 at the C-terminus.		
Source:	Baculovirus-Insect cells		
Molecular Weight:	The recombinant human TNFRSF11B consists of 618 amino acids and predicts a molecular mass of 70.3 kDa.		
Formulation:	Lyophilized from sterile pH 7.0, 100 mM Glycine, 10 mM NaCl. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization.		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	< 1.0 EU/µg as determined by LAL method.		
Predicted N terminal:	Glu 22		
Reconstitution:	<p>Centrifuge vial prior to opening. Add sterile distilled water to a concentration of 0.1 mg/mL and gently pipette the solution up and down the sides of the vial.</p> <p>DO NOT VORTEX. Allow several minutes for complete reconstitution.</p>		
Storage & Stability:	<p>Stable for up to 1 year from date of receipt at -20°C to -80°C</p> <p>After reconstitution, store working aliquots at -20°C to -80°C.</p> <p>Avoid repeated freeze-thaw cycles.</p>		



SDS-PAGE



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com