

## SERPINE1

### Recombinant Cynomolgus Monkey PAI-1 Elastase Cleaved

<b>Catalog No.</b>	CRK004A CRK004B	<b>Quantity:</b>	500 µg 1.0 mg
<b>Description:</b>	Recombinant cynomolgus monkey wild type PAI-1 is produced by cleavage at P3-P4 residues via immobilized elastase.		
<b>Concentration:</b>	1.0 mg/ml		
<b>Gene ID:</b>	713883		
<b>Source:</b>	<i>E. coli</i>		
<b>Molecular Weight:</b>	39 kDa		
<b>Formulation:</b>	Frozen liquid in 0.05 M sodium phosphate + 0.1 M NaCl + 1 mM EDTA, pH 6.6.		
<b>Biological Activity:</b>	<1% active		
<b>Applications:</b>	The elastase cleaved monkey PAI-1 is useful for various control experiments requiring a non-reactive, reactive loop inserted species of PAI-1.		
<b>Storage &amp; Stability:</b>	Store at -80°C. Stable for 3 years from delivery. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

