

LLO

Recombinant Human Listeriolysin-O

Catalog No.	CRL121A	Quantity:	10 µg
	CRL121B		50 µg
	CRL121C		100 µg

Alternate Names: Listeriolysin-O, LLO, hlyA.

Description: Listeriolysin O (aka LLO) is a hemolysin produced by *Listeria monocytogenes* bacteria, the pathogen responsible for causing listeriosis. The toxin may be regarded as a virulence factor, since it is crucial for the virulence of *L. monocytogenes*. LLO is a single polypeptide protein encoded by the hlyA gene and composed of 529 residues. LLO is a thiol-activated cholesterol-dependent pore forming toxin protein; therefore, it is activated by reducing agents and inhibited by oxidizing agents. Still, LLO differs from other thiol-activated toxins, as its cytolytic activity is maximized at a pH of 5.5. Inside the acidic phagosomes (average pH ~ 5.9) of cells that have phagocytosed *L. monocytogenes*, LLO is selectively activated by maximizing activity at a pH of 5.5. Following the phagosome lysis by LLO, the bacterium breaks out into the cytosol, where it is able to grow intracellularly, and the toxin has reduced activity in the more basic cytosol. Thus, LLO permits *L. monocytogenes* to break out from the phagosomes into the cytosol without harming the plasma membrane of the infected cell, which allows the bacteria to live intracellularly, where they are sheltered from extracellular immune system factors such as the complement system and antibodies. LLO also brings about dephosphorylation of histone H3 and deacetylation of histone H4 in the early phases of infection, before entry of *L. monocytogenes* into the host cell. The pore-forming activity is not implicated in causing the histone modifications. The modifications of the histones affect the down regulation of genes encoding proteins involved in the inflammatory response. Therefore, LLO may be significant in subverting the host immune response to *L. monocytogenes*. At its NH₂-terminus it possesses a 25 residues long typical signal sequence excised during the secretion process. Moreover, in its NH₂-terminus there is also a 19 amino acids PEST-like sequence that may target this toxin for degradation. The PEST-like sequence found in LLO and is considered crucial for virulence, given that mutants lacking the sequence lysed the host cell. Nevertheless, contrary to PEST's supposed role in protein degradation, evidence implies that the PEST-like sequence may control LLO production in the cytosol rather than increase degradation of LLO. This hemolysin protein is provided in its mature form that is constituted by 503 amino acids, from 26 to 529 of the registered sequence (accession number: AAF64524).

Physical Appearance: Sterile Filtered clear solution.

Protein Accession No: AAF64524

Source: *E. coli*

Molecular Weight: 58 kDa

Formulation:



The protein contains 50 mM NaH₂PO₄ + 1 mM EDTA + 2.7 mM KCl + 1 mM DTT + 5% glycerol and 0.5 M NaCl.

- Activity:** LLO has a hemolytic activity 1.25×10^6 HU/mg.
2 mM DTT is used to reactivate the toxin.
- Amino Acid Sequence:** This hemolysin protein is provided in its mature form that is constituted by 503 amino acids, from 26 to 529 of the registered sequence (accession number: AAF64524).
- Applications:**
- Cytosolic delivery of molecules, peptides, oligonucleotides and plasmid DNA.
 - Production of specific monoclonal antibodies.
 - Detection of anti-listeriolysin O antibodies.
- Storage & Stability:** LLO although stable at room temperature for 1 week, should be stored desiccated below -18°C.
Please prevent freeze-thaw cycles.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

