

Animal, Bacterial & Viral Free - Low Endotoxin - Ultra Pure - High BioActivity
FGF2

Recombinant Human Fibroblast Growth Factor 2, FGF2, Endotoxin Free

Catalog No.	CRO500A CRO500B CRO500C	Quantity:	10 µg 50 µg 100 µg
Alternate Names:	FGF basic (146 aa), FGF-2, HBGF-2, Prostatropin		
Description:	Fibroblast Growth Factor 2 is a member of the FGF family heparin-binding multifunctional proteins. These proteins play a role in wound healing, angiogenesis, and embryonic development, during pre- and post-natal growth and regeneration of tissues, by promoting cellular proliferation and differentiation. FGF 2 is recognized as a critical component of human embryonic stem cell culture medium due to its vital role for cells to remain in an undifferentiated state. Recombinant human FGF 2 consists of 146 aa residues and a 10 aa His-tag at the N-terminus.		
UniProtKB:	P09038		
Gene ID:	2247		
Source:	<i>Hordeum vulgare</i> (barley grain). Barley grain's proteolytic activity is almost 50 times less than <i>E. coli</i> or mammalian cells. Barley seed has no human or animal viral contaminants, which is ideal for stem cell culture and <i>in vitro</i> and <i>in vivo</i> studies.		
Molecular Weight:	17.8 kD		
Formulation:	Lyophilized from a 0.2 µm sterile filtered solution of PBS, pH 7.2.		
Purity:	>95% by SDS-PAGE gel analysis.		
Endotoxin Level:	< 0.005 ng per µg of product (0.05 EU/µg) as measured by turbidimetric kinetic assay.		
Biological Activity:	ED ₅₀ < 0.4 ng/ml. Bioactivity was determined by dose dependent effect on proliferation of 3T3 cells.		
Specific Activity:	> 2.5 x 10 ⁶ U/mg		
Reconstitution:	Centrifuge vial prior to opening. It is recommended to reconstitute the lyophilized protein in sterile water to a concentration of no less than 100 µg/ml. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please note that the addition of any carrier protein into this product may introduce endotoxin. Depending upon the particular application employed, this may be undesirable.		
Storage & Stability:	The lyophilized protein, though stable at room temperature for two weeks, is best stored at -20°C. Reconstituted protein should be used immediately or stored in working aliquots at -20°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

