

Animal, Bacterial & Viral Free - Low Endotoxin - Ultra Pure - High BioActivity
BMP2

Recombinant Human Bone Morphogenetic Protein 2, BMP2, Endotoxin Free

Catalog No.	CRO506A CRO506B CRO506C CRO506D	Quantity:	10 µg 50 µg 1 mg 100 µg
Alternate Names:	BMP-2A, BMP2A		
Description:	Human Bone Morphogenetic Protein 2 (BMP-2) is a member of the TGF-beta superfamily of cytokines that effect bone and cartilage formation. Human BMP-2 is a disulphide linked homodimeric protein. BMP-2 plays an important role in various processes during embryonic development, including induction of bone and cartilage growth and regeneration, cardiac cell differentiation and epithelial to mesenchymal transition. BMP2 promotes osteogenic and chondrogenic differentiation from mesenchymal cells, stem cells, and vascular smooth muscle cells. BMP-2 has great potentials for regenerative and medical therapeutic applications.		
Gene ID:	650		
UniProtKB:	P12643		
Source:	Recombinant Human BMP-2 is produced in <i>Nicotiana tabacum</i> that is void of any human or animal viral contaminants that could jeopardize stem cell culture.		
Molecular Weight:	Recombinant human BMP-2 produced in transgenic tobacco plants is a homodimeric polypeptide of 114 aa with a MW of 26 kDa		
Formulation:	Lyophilized from a sterile solution of 25 mM MES, pH 6.0		
Purity:	> 95% by SDS-PAGE gel analysis.		
Endotoxin Level:	Endotoxin level is less than 0.005 ng per µg (0.05 EU/µg) as measured by kinetic LAL assay.		
Biological Activity:	ED ₅₀ < 100 ng/ml, determined by the dose dependent induction of alkaline phosphatase production by ATDC-5 mouse chondrogenic cells.		
Specific Activity:	> 1.0 x 10 ⁴ U/mg		
Reconstitution:	Centrifuge vial prior to opening. Reconstitute the lyophilized protein in sterile water to a concentration of no less than 50 µg/ml. Due to the protein nature, dimers and multimers may be observed. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please note that the addition of any carrier protein into this product may introduce endotoxin. Depending upon the particular application employed, this may be undesirable.		
Storage & Stability:	The lyophilized protein, though stable at room temperature for few weeks, is best stored at -20°C. Reconstituted protein should be used immediately or stored in working aliquots at -20°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

