

Recombinant Human Procalcitonin/CALCA His

Catalog No.	CRP143A	Quantity:	2 µg
	CRP143B		10 µg
	CRP143C		1.0 mg

Alternate Names: PCT

Description: Recombinant Human Procalcitonin is a single, non-glycosylated, polypeptide chain containing 114 amino acids fragment (3-116) and an amino-terminal hexahistidine tag.

Source: *E. coli*

Molecular Weight: 17.13 kDa

Formulation: Sterile filtered liquid in PBS + 50% glycerol

Purity: > 95.0% as determined by RP-HPLC and SDS-PAGE analyses

Endotoxin Level: < 0.1 ng/µg of PCT

Storage & Stability: Stable for 4 weeks at 2-4°C or in working aliquots at -20°C for longer storage. **Avoid repeated freeze-thaw cycles.**

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

