

Plg Native Rat Plasminogen

Catalog No.	CRP260B CRP260C	Quantity:	1 mg 10 mg
Description:	Native rat Plasminogen prepared from fresh rat plasma by immobilized lysine chromatography.		
UniProt ID:	Q01177		
Gene ID:	85253		
Source:	Rat plasma		
Molecular Weight:	92.0 kDa		
Formulation:	Liquid in 0.1 M HEPES, 0.1 M NaCl, pH 7.4		
Purity:	>95% by SDS-PAGE analysis		
Concentration:	≥ 1.0 mg/ml, lot specific $E^{0.1\%}_{280\text{nm}} = 1.69$		
Storage & Stability:	Store at -80°C for up to 1 year. Upon initial thaw, prepare working aliquots and store at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.