

PLAT

Recombinant Human tPA FITC

Catalog No.	CRT172B CRT172C	Quantity:	100 µg 1.0 mg
Alternate Names:	Tissue Plasminogen Activator, T-Plasminogen Activator, PLAT, T-PA		
Description:	Recombinant Human Tissue Plasminogen Activator FITC labeled, > 85% single chain. 100% complex formation with PAI-1.		
UniProt ID:	P00750		
Gene ID:	5327		
Source:	CHO cells		
Molecular Weight:	70.0 kDa		
Formulation:	Frozen liquid in 0.4 M HEPES, 0.1 M NaCl, pH 7.4		
Purity:	>95% by SDS-PAGE		
Concentration:	~ 2 mg/ml, lot specific		
F/P Ratio:	Lot specific		
Extinction Coefficient:	1.9		
Storage & Stability:	Store in working aliquots, protected from light, at -80°C for up to 1 year. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

