

## TNFSF10

### Recombinant Human soluble TRAIL/APO2L

<b>Catalog No.</b>	CRT500A CRT500B CRT500C	<b>Quantity:</b>	10 µg 50 µg 1.0 mg
<b>Alternate Names:</b>	TNF-related apoptosis-inducing ligand, Apo-2 ligand, Tumor necrosis factor ligand superfamily member 10, TNFSF10, CD253		
<b>Description:</b>	TNF-Related Apoptosis-Inducing Ligand (TRAIL) is a cytokine produced and secreted by most normal tissue cells. It is shown to be a cytotoxic protein that induces apoptosis in tumor cells through activation of the death receptors, DR4 and DR5. TRAIL and its receptors have been used as the targets of several anti-cancer therapeutics since the mid-1990s, such as Mapatumumab. However, these have not shown significant survival benefit. TRAIL has also been implicated as a pathogenic or protective factor in various pulmonary diseases, particularly pulmonary arterial hypertension. Human TRAIL is active on mouse cells.		
<b>Gene ID:</b>	8743		
<b>UniProt ID:</b>	P50591		
<b>Source:</b>	<i>E. coli</i>		
<b>Molecular Weight:</b>	19.5 kDa (168 aa)		
<b>Formulation:</b>	Lyophilized from a sterile filtered solution of 10 mM sodium phosphate, 50 mM NaCl, pH 7.5		
<b>Purity:</b>	≥ 90.0% by reducing and non-reducing SDS-PAGE		
<b>Endotoxin Level:</b>	≤1 EU/µg by kinetic LAL analysis.		
<b>Biological Activity:</b>	Typical ED <sub>50</sub> < 15 ng/ml, determined by the ability to induce apoptotic cell death in TRAIL-sensitive RPMI-8226 cells.		
<b>Amino Acid Sequence:</b>	MRERGPQRVA AHITGTRGRS NTLSSPNSKN EKALGRKINS WESSRSGHSF LSNLHLRNGE LVIHEKGFYY IYSQTYFRFQ EEIKENTKND KQMVQYIYKY TSYPDPILLM KSARNSCWSK DAEYGLYSIY QGGIFELKEN DRIFVSVTNE HLIDMDHEAS FFGAFLVG		
<b>Reconstitution:</b>	<b>Centrifuge vial prior to opening.</b> When reconstituting the product, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution. <b>Do not vortex.</b> It is recommended to reconstitute the lyophilized product with sterile water at a concentration of 0.1 mg/ml. After complete solubilization of the protein, it can be further diluted to other aqueous solutions.		



**Storage & Stability:**

Lyophilized product is stable at room temperature for shipping purposes. Upon receipt, store at -20°C to -80°C for up to 1 year.

Upon reconstitution, the preparation is stable for up to one month at 2-8°C. For long term storage, prepare working aliquots and store at -20 to -80°C. For maximal stability, dilute to working aliquots in a 0.1% BSA solution.

**Avoid repeated freeze-thaw cycles.**


**Human TRAIL Gel**

Figure: 1 ug run under (-) non-reducing conditions and (+) reducing conditions in a 4-20% Tris-Glycine gel, stained with Coomassie Blue. Human TRAIL is predicted to have a MW of 19.5 kDa.

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

