

## TSLP

### Recombinant Human Thymic Stromal Lymphopoietin

<b>Catalog No.</b>	CRT702A CRT702B CRT702C	<b>Quantity:</b>	2 µg 10 µg 1.0 mg
<b>Alternate Names:</b>	TSLP		
<b>Description:</b>	Thymic Stromal Lymphopoietin (TSLP) is a hemopoietic protein that is expressed in the heart, liver and prostate. TSLP overlaps biological activities with IL-7 and binds with the heterodimeric receptor complex consisting of the IL-7R alpha chain (IL-7R alpha) and the TSLP-specific chain (TSLPR). Like IL-7, TSLP induces phosphorylation of STAT3 and STAT5, but uses kinases other than the JAKs for activation. TSLP prohibited apoptosis and stimulated growth of the human acute myeloid leukemia (AML)-derived cell line MUTZ3. It induces the release of T cell-attracting chemokines TARC and MDC from monocytes and activates CD11c(+) dendritic cells (DCs). TSLP activated DCs primed naive T cells to produce the proallergic cytokines (IL-4, IL-5, IL-13, TNFalpha) while down-regulating IL-10 and IFN-gamma suggesting a role in initiating allergic inflammation.		
<b>Physical Appearance:</b>	Sterile Filtered White lyophilized (freeze-dried) powder.		
<b>Gene ID:</b>	85480		
<b>Source:</b>	<i>E. coli</i>		
<b>Molecular Weight:</b>	~ 15 kDa		
<b>Formulation:</b>	Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM PB, pH 7.4 + 150 mM NaCl.		
<b>Purity:</b>	> 98% as determined by SDS-PAGE and HPLC analyses		
<b>Endotoxin Level:</b>	Less than 1EU/µg of rHuTSLP as determined by LAL method.		
<b>Biological Activity:</b>	Fully biologically active when compared to standard. The ED <sub>50</sub> determined by a cell proliferation assay using human IL-7R alpha and human TSLP R co-transfected murine BaF3 pro-B cells is less than 0.3 ng/ml.		
<b>Amino Acid Sequence:</b>	MYDFTNCDFE KIKAAYLSTI SKDLITYMSG TKSTEFNNTV SCSNRPHCLT EIQLTFNPT AGCASLAKEM FAMKTKAALA IWCPGYSETQ INATQAMKKR RKRKVTTNKC LEQVSQLQGL WRRFNRPLLK QQ		
<b>Reconstitution:</b>	<b>Centrifuge vial prior to opening.</b> Add sterile distilled water or aqueous buffer to a concentration of 0.1-1.0 mg/mL. Further dilutions should be made in appropriate buffered solutions.		
<b>Storage &amp; Stability:</b>	This lyophilized preparation is stable at 2-4°C, but should be kept desiccated at -20°C for long term storage. Upon reconstitution, the preparation is stable for up to one week at 2-4°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -80°C. <b>Avoid repeated freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.