

PLAU

Recombinant Human LMW Urokinase HRP

Catalog No.	CRU116B CRU116C	Quantity:	100 µg 1.0 mg
Alternate Names:	Urokinase Plasminogen Activator, PLAU, ATF, UPA, URK, u-PA		
Description:	Horseradish peroxidase conjugated LMW urokinase purified from auto-catalytic digestion of HMW human urokinase. LMW uPA contains the proteinase domain and is separated from the amino terminal fragment by cleavage at the Lys135-Lys136 bond.		
Concentration:	2.5 mg/ml		
Gene ID:	5328		
Source:	Insect cell culture		
Molecular Weight:	33 kDa		
Formulation:	Liquid in 0.05 M sodium acetate + 0.1 M NaCl + 1 mM EDTA + 50% glycerol, pH 5.0		
Storage & Stability:	Store in working aliquots at -20°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

