

Recombinant Mouse HMW Urokinase FITC

Catalog No.	CRU120A CRU120B	Quantity:	50 µg 100 µg
Alternate Names:	Urokinase Plasminogen Activator, PLAU, ATF, UPA, URK, u-PA		
Description:	Recombinant Mouse High Molecular Weight Urokinase FITC labeled		
Concentration:	2.0 mg/ml		
Gene ID:	18792		
Protein Accession No:	NP_032899		
Source:	Insect cells		
Molecular Weight:	44.0 kDa		
Formulation:	Liquid in 0.05 M sodium acetate + 0.1 M NaCl + 1 mM EDTA, pH 5.0		
Storage & Stability:	Store in working aliquots, protected from light, at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

