

VHL

Recombinant Human Von-Hippel-Lindau Tumor Suppressor Protein His

Catalog No.	CRV108A CRV108B CRV108C	Quantity:	10 µg 50 µg 1.0 mg
Alternate Names:	Von Hippel-Lindau disease tumor suppressor, pVHL, Protein G7, VHL, RCA1, VHL1, HRCA1.		
Description:	<p>Von Hippel-Lindau disease is a dominant inherited syndrome characterized by the predisposition to develop various kinds of benign and malignant tumors, including clear cell renal carcinomas, pheochromocytomas and hemangioblastomas of the central nervous system and retina. VHL syndrome is caused by germline mutation in the VHL tumor suppressor, and VHL tumors are associated with loss or mutation of the remaining wild-type allele. VHL has two domains: a roughly 100-residue NH₂-terminal domain rich in b-domain and a smaller alpha-helical domain (a-domain), held together by two linkers and a polar interface. VHL protein is also involved in the degradation of hypoxia-inducible factor (HIF).</p> <p>Recombinant Human Von Hippel-Lindau Protein b-domain is a single, non-glycosylated polypeptide chain containing aa 1-154 fused to a 20 aa His-Tag at the N-terminus for a total of 174 aa.</p>		
Physical Appearance:	Sterile filtered colorless solution.		
GenelD:	7428		
Source:	<i>E. coli</i>		
Molecular Mass:	19.2 kDa		
Formulation:	The Von Hippel-Lindau Protein contains PBS + pH-7.4 + 2 mM EDTA + 1 mM DTT.		
Purity:	Greater than 95.0% as determined by SDS-PAGE.		
Purification:	Purified by proprietary chromatography techniques.		
Amino Acid Sequence:	MGSSHHHHHH SSGLVPRGSH MPRRAENWDE AEVGAE EAGV EEYGPEEDGG EESGAEESGPEESGPEELGA EEEMEAGRPR PVLRSVNSRE PSQVIFCNRS PRVLPVWLN FDGEPQPYPT LPPGTGRRIH SYRGHLWLFR DAGTHDGLLV NQTELFVPSL NVDGQPIFAN ITLP.		
Storage & Stability:	Store at 4°C if entire vial will be used within 2-4 weeks. Store frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

