

VTN

Native Human Monomeric Vitronectin

Catalog No.	CRV123B CRV123C	Quantity:	500 µg 1 mg
Alternate Names:	V75, VN, VNT		
Gene ID:	7448		
Protein Accession No:	NP_000629		
Description:	Native human monomeric Vitronectin prepared from fresh human plasma using non-denaturing chromatography		
Concentration:	1.0 mg/ml		
Source:	Human plasma		
Molecular Weight:	75.0 kDa		
Formulation:	Liquid in 0.05 M sodium phosphate + 0.1 M NaCl, pH 7.4		
Storage & Stability:	Store in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

