

## Vtn

### Native Mouse Multimeric Vitronectin

<b>Catalog No.</b>	CRV128B CRV128C	<b>Quantity:</b>	100 µg 1 mg
<b>Alternate Names:</b>	V75, VN, VTN, S-protein, Serum-spreading factor		
<b>Description:</b>	Native mouse multimeric Vitronectin prepared from fresh mouse plasma using urea as a denaturant.		
<b>UniProt ID:</b>	P29788		
<b>Gene ID:</b>	22370		
<b>Source:</b>	Mouse plasma		
<b>Molecular Weight:</b>	75 kDa monomer		
<b>Formulation:</b>	0.05 M sodium phosphate, 0.1 M NaCl, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE analysis		
<b>Concentration:</b>	≥ 1.0 mg/ml, lot specific		
<b>Extinction Coefficient:</b>	$E^{0.1\%}_{280\text{nm}} = 1.38$		
<b>Storage &amp; Stability:</b>	Store at -80°C for up to 1 year. Upon initial thaw, prepare working aliquots and store at -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

