

Serpinc1

Rabbit Anti Mouse Antithrombin III pAb

Catalog No.	CS190A CS190B	Quantity:	1.0 ml 10 ml
Alternate Names:	At3, At-3, ATIII, AI114908		
Description:	Polyclonal antiserum (host rabbit) to mouse antithrombin.		
Concentration:	Whole Antiserum		
Gene ID:	11905		
Purity:	Serum		
Specificity:	Recognizes Antithrombin III		
Host:	Rabbit		
Immunogen:	Full length native protein (purified) (Mouse).		
Isotype:	IgG		
Formulation:	Liquid in Whole Antiserum		
Cross-Reactivity:	Mouse		
Applications:	ELISA, WB		
Storage & Stability:	When stored at -80°C, product is stable for 5 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

