

Ren1

Rabbit Anti-Mouse Prorenin/Renin IgG fraction FITC pAb

Catalog No.	CS206A CS206B	Quantity:	1.0 mg 10 mg
Alternate Names:	Ren, Rnr, Rn-1, Ren-1, Ren-A, Ren1c, Ren1d, D19352		
Description:	FITC labeled polyclonal antibody (host rabbit) to mouse prorenin. IgG fraction purified by immobilized Protein A.		
Concentration:	10.0 mg/ml		
Gene ID:	19701		
Purity:	IgG fraction		
Specificity:	Recognizes Renin/Prorenin		
Host:	Rabbit		
Immunogen:	Recombinant full length mouse protein fused to a His tag at C-terminal end (Glycosylated).		
Isotype:	IgG		
Conjugate:	FITC labeled		
Formulation:	Liquid in 0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6		
Purification:	IgG fraction purified by immobilized Protein A.		
Cross-Reactivity:	Mouse		
Applications:	ELISA, WB		
Storage & Stability:	When stored at -80°C, product is stable for 5 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

